

# UVICO – Continuous Light Source



- ▶ life time 2000 hours
- ▶ 120 W UVICO 2000 lamp, pre-aligned, exchange without adjustment
- ▶ fibre optically coupled, broad range of optical fibres available
- ▶ built-in fast shutter, manual and remote (TTL) control
- ▶ shutter opening time 30ms -  $\infty$
- ▶ filter slider with 3 positions
- ▶ input 110 - 240 V ac, 50/60 Hz
- ▶ options: remote control via Ethernet featuring
  - ◆ ON / OFF ◆ shutter ◆ attenuator for 1% to 100% ◆ filter changer

